


## 2680A-OPC | OPC software for Fluke 2680 Series



### Product overview: 2680A-OPC | OPC software for Fluke 2680 Series

OPC is a standard interface for passing data and communicating information between programs and devices in the Windows environment. When using the 2680A-OPC software, client programs such as National Instruments LabVIEW™, Canary Labs Trend Link, Wonderware® software and many others may request real-time or historical data from the Fluke DAQ database.

#### Requirements

700 MHz processor or greater

Windows 2000 or later

Free disk space – 400 MB

## Ordering information



### **2680A-OPC**

OPC softWare for Fluke 2680 Series

---



**Fluke.** *Keeping your world up and running.®*

**Fluke Corporation**  
PO Box 9090, Everett, WA 98206 U.S.A.

**For more information call:**  
In the U.S.A. (800) 443-5853  
In Canada (800) 36-FLUKE  
From other countries +1 (425) 446-5500  
[www.fluke.com](http://www.fluke.com)

©2025 Fluke Corporation.  
Specifications subject to change without notice.  
04/2025

**Modification of this document is not permitted  
without written permission from Fluke  
Corporation.**